


Praesto[®]

Jetted 35 µm

Ion Exchange
Chromatography Resins



Your solutions company


Purolite Life Sciences focuses on any applications involving interactions with people, bringing Purolite's innovative thinking and distinguished history of resin technology expertise to the global Life Sciences marketplace.

APIs, enzyme carriers, immobilized enzymes, and agarose or synthetic chromatography resins for purification and separation, to support research and development and production-scale applications in pharmaceuticals, protein purification, food processing, bioprocessing and fine chemical markets.

**“We provide solutions
for our customers’
most critical questions.”**

The background of the slide is a photograph of a clear blue sky with two large, fluffy white clouds. One cloud is in the bottom left corner, and the other is in the bottom right corner, both appearing to rise towards the center of the frame.

**Your trusted partners
for resin technology
solutions.**



Jetting - A patented,
innovative manufacturing
process for **superior
resin bead uniformity,
performance and column
packing consistency.**

Introduction

Purolite® Life Sciences has revolutionized the future of Protein A chromatography with Praesto® Jetted resins. It is the only bioprocess-scale agarose resin available with uniform particle sizes for superior performance.

Manufactured with highly crosslinked agarose, Praesto SP 35 µm (cation) and Praesto Q 35 µm (anion) provide high resolution polishing resins for recombinant proteins and other biomolecules.

Key features



Uniform particle sizes



Increased dynamic binding capacity



Chemically stable



Supplied in bulk or in pre-packed, pre-qualified OPUS® columns from Repligen Corporation



Non-solvent, green manufacturing

Praesto® Jetting Advantages

Advantages
of Jetting

The diagram consists of seven circles arranged in two rows. The top row has three circles, and the bottom row has four circles. The first circle in the top row is shaded gray and contains the text 'Advantages of Jetting'. The other six circles are white with gray outlines and contain specific advantages in blue text. The background is a light gray grid of small dots.

**More consistent
packing characteristics**


**Better
kinetics**

**Improved
resolution**

**Higher dynamic
binding capacities**

**Reduced buffer
consumption**

**Longer resin
lifetime**




**“Purolite Life Sciences
is the only provider
in the world capable
of supplying these
enhanced, uniform
beads in process-scale
volumes.”**

Background

Praesto® ion exchange resins are designed for biomolecule purification, including proteins, peptides and oligonucleotides.

They are manufactured using 'Jetting' - a new, patented manufacturing method that produces agarose beads with a very narrow particle size distribution. Purolite is the first agarose resin provider to produce process-scale volumes of these uniform particle size beads.



Praesto® Jetted resins demonstrate superior performance characteristics including:

- Increased pressure/flow properties
- Improved resolution from uniform particles
- Improved packing reproducibility and stability
- Higher dynamic binding capacity
- Green manufacturing

Figure 1: Microscope image - Praesto® Jetted agarose resins

A hand holding a glowing lightbulb against a bokeh background. The background is dark with many out-of-focus light spots in shades of green and yellow. The lightbulb is held in the center-right of the frame, and its glow illuminates the hand and the surrounding bokeh.

Jetting removes the need for the extensive sieving process found in traditional batch wise emulsification technology.


This generates considerably less waste and results in shorter lead times.

The Jetting Process


Praesto® Jetted products are manufactured using a patented, continuous technology called 'Jetting'.

This process involves an agarose solution pumped through a specially-engineered steel membrane or can into a mineral oil, forming a suspension with a narrower particle size distribution (UC <1.3) than traditional resins (UC <1.9).

A high yield, economical process, it provides Praesto® Jetted products with unique performance characteristics."



The Jetting Process continued...



S. R. Kosvintsev, Pat. 20140264984, 2013, Method of producing uniform polymer beads of various sizes. US Pat. 9028730 B2, 2015

Environmental Benefits


Jetted uniform beads also have several important benefits to the environment over standard agarose resin beads. Jetting technology is a continuous manufacturing process that produces beads in a more efficient way, with very high yields, dramatically lowering lead times and improving the overall supply chain.

Since the beads are uniform in size, it removes the need for extensive sieving, thus generating considerably less waste. Another environmentally friendly advantage with this new technology, is that it eliminates the need for high levels of organic solvents compared to standard agarose resin beads, like toluene commonly used in traditional manufacturing and is a major advancement in 'green manufacturing'



ENVIRONMENTALLY
FRIENDLY

JETTING
TECHNOLOGY



HIGH
YIELDS



The Praesto® Range

The Praesto® range offers a selection of modern, high-flow Affinity and Ion Exchange agarose resins, delivering exceptional results from Protein A to high-resolution polishing steps. The range also includes a full selection of Praesto® Pure base matrices, and pre-activated resins in a variety of source chemistries.

All Praesto® products provide an advanced, high-flow, highly cross-linked agarose base matrix. The entire range benefits from excellent pressure/flow characteristics and stability for optimal recovery of active proteins.



Praesto® products are also available pre-packed in OPUS® columns from Repligen Corporation.

The perfect choice for screening, sample preparation, and process validation through scale-up and commercial manufacturing.



Ligand Structure

The strong ion exchange ligand groups of Praesto SP (a) and Praesto Q (b) are well established in large scale purification.




Figure 2. Ion structure of Praesto SP and Praesto Q

Praesto® Q: Typical Physical & Chemical Characteristics

Application

mAb purification

Polymer Structure

Cross-linked agarose,
uniform size beads

Appearance

Spherical, uniform size beads
supplied in 20% ethanol slurry.

Functional Group

$\text{CH}_2\text{N}^+(\text{CH}_3)_3$

Binding Capacity

> 80 mg/ml BSA at 6 mins resi-
dence time

Average Particle Size

35 μm

Ionic Capacity, mmol/mL Resin

0.14 - 0.18

pH Stability, (Long-term)

3-13

pH Stability (Short-term)

2-14

Recommended Storage

20% ethanol, 4-30°C

Flow velocity

>150

1 d50v is the median particle size of the cumulative volume distribution.

Praesto® SP: Typical Physical & Chemical Characteristics

Application

mAb purification

Binding Capacity

> 90 mg/ml IgG at 6 mins
residence time

pH Stability (Short-term)

3-14

Polymer Structure

Cross-linked agarose,
uniform size beads

Average Particle Size

35 µm

Recommended Storage

20% ethanol, 4-30°C

Appearance

Spherical, uniform size beads
supplied in 20% ethanol slurry.

Ionic Capacity, mmol/mL Resin

0.11 - 0.16

Flow velocity

>150

Functional Group



CH₂CH₂CH₂SO₃⁻

pH Stability, (Long-term)

4-13

1 d50v is the median particle size of the cumulative volume distribution.

Particle Size Distribution



The particle size distribution of jetted agarose is significantly improved compared to batch emulsification.

A 35 μm jetted = 95% within 25 – 50 μm with a uniformity coefficient* of <1.3.

A 45 μm batch emulsified = 95% within 25 - 100 μm with a uniformity coefficient of 1.9.

Figure 3. Particle Size Distribution

*Uniformity coefficient (UC) is defined as the $D_{60}(V)$ divided by the $D_{10}(V)$.

Pressure Flow

With new crosslinking methods, combined with jetting technology, Purolite Life Sciences have been able to develop 35 µm small particles for improved resolution for high polishing applications that are still suitable for large scale bioprocessing. For instance, a 45 µm batch emulsified bead has a lower pressure flow rate capability compared to a 45 µm jetted bead. When packed at 4 bar and ran at 3 bar, in a HiScale 26/40 (with a 2.6 cm diameter and a 20 cm BH), the 45 µm batch emulsified resin flow velocity reached 220 cm/h compared to the same size jetted resin (45 µm) which reached a maximum flow velocity of 282 cm/h.



D ₅₀	IgG DBC (10%) 2.4 min	IgG DBC (10%) 6 min	Flow Velocity (cm/h) 20 cm bed height
Batch Emulsification	68	97	220
45 (Jetted)	80	108	282


Application Example

To demonstrate the high resolution of jetted agarose base beads over batch emulsified resins, we assessed the separation profiles of Praesto[®] cation exchangers to reduce the acidic species using a customer-supplied biosimilar MAb.

In this example, the goal was to ensure the acidic species of the biosimilar matched that of the originator. Purification with Praesto[®] AP protein A affinity resin was followed by an anion multimodal in flow mode as the first polishing step. The material was then further purified in a bind and elution mode using two different CIEX resins.

- 1) SP 45 μm - a traditional batch emulsified agarose resin (95% within 25-100 μm)
- 2) Jetted SP 35 μm , a uniform particle size agarose resin (95% within 25-50 μm , UC <1.3).

The separation profiles were assessed using a standard conductivity gradient elution.



Comparing SP 45 (batch emulsified) with Jetted SP 35, a far more pronounced shoulder of the acidic species can be seen on the chromatograms.




Figure 4. Comparing SP 35 and SP 45



A photograph of two men in a warehouse setting. The man on the left is wearing a white hard hat, a light blue button-down shirt, and a patterned tie. He is holding a tablet and looking at it. The man on the right is wearing a yellow hard hat, a dark blue polo shirt, and a high-visibility yellow safety vest. He is gesturing with his hands while looking at the tablet. In the background, there are stacks of blue and green industrial drums on wooden pallets. The text "Stability from a robust supply chain" is overlaid in the top right corner in a white, sans-serif font.

Stability from
a robust supply
chain

Complete regulatory support

For over 35 years, Purolite® has supplied specialty resin technology to industries within complex regulatory environments including biotechnology, pharmaceutical, food, fine chemical and electric power generation.

The regulatory environment is ever changing, driven by increasing regulatory requirements, increasing development costs and times, and market pressures impacting pharma and food industries.

For Life Sciences products, Regulatory Support Files (RSF) are available. Regulatory Support Files provide direct and detailed information on performance, stability, extractable compounds, and analytical methods for each resin.



“Regulatory expertise throughout the product life-cycle is essential to identify options for product development, optimize ‘speed to market’ and produce a product that meets customer needs.”

Purolite® implements control documentation and processes at every level to ensure regulatory support to customers using our products.

Purolite® complies with required national and international regulations, as well as many voluntary specialty certifications.

These include:

- GMO/TSE/BSE free
- ISO 9001:2015 quality system specifications
- ISO 14001:2015 Environmental Management System requirements

We also hold Drug Master Files with the US FDA, Japan, Canada and EU



**Security
of supply**

**Lot-to-lot
consistency**

**Consistent operating
practices across sites**

**Internal
audits**

**Continual quality
training**

**Compliant
containers**

**Certificates of
analysis (CoA)**

**Raw material
supplier program**

Safety data

Purolite Life Sciences maintains Material Safety Data Sheets (MSDS for the U.S. and ERSDS for Europe) on each of its bulk resins. These data sheets contain relevant information that you may need to protect your customers and employees against any known health or safety hazards associated with our products.

Purolite Life Sciences supplies copies of our Material Safety Data Sheets with all bulk resins. These describe precautions to be taken in the storage and handling of our products and in the maintenance of the health and safety of persons exposed to our products, the public and the environment with respect to our products.



Placing your order



How to order

To place your order simply contact us via email or telephone using the information on the next page, and quote your order number from the table below. Praesto® Jetted A50 can be supplied loose in the following pack sizes; 25ml, 100ml, 500ml, 1L+ on request.

It can also be supplied pre-packed into OPUS® columns by Repligen®. Available in RoboColumns with 200 µl column volumes, and MiniChrom available in 1 and 5 ml column volumes.

For scale up/validation, pilot manufacturing and cGMP manufacturing operations we have an agreement with Repligen® to pack Praesto® resins in 0.5 cm diameter to 80 cm diameter OPUS® columns with flexible bed heights.

If you would like to discuss how Praesto® Jetted SP or Q can benefit your purification process, we have dedicated experts on-hand across the globe to provide knowledgeable, same-day technical assistance.

How to order - Jetted SP 35

BULK RESIN

PACK SIZE

ORDER NUMBER

Praesto® Jetted SP35	25 ml	PR00432-166
Praesto® Jetted SP35	100 ml	PR00432-164
Praesto® Jetted SP35	500 ml	PR00432-165
Praesto® Jetted SP35	1 L	PR00432-310
Praesto® Jetted SP35	5 L	PR00432-311
Praesto® Jetted SP35	10 L	PR00432-312

PRODUCT

COLUMN SIZE

ORDER NUMBER

Praesto® Jetted SP35 HT	1 ml Column	PR00432-275*
Praesto® Jetted SP35 HT	5 ml Column	PR00432-276*

** 1 ml HT columns packed with Ion Exchange resins available in packs of 5 only*

How to order - Jetted Q 35

BULK RESIN

PACK SIZE

ORDER NUMBER

Praesto® Jetted Q35	25 ml	PR00436-166
Praesto® Jetted Q35	100 ml	PR00436-164
Praesto® Jetted Q35	500 ml	PR00436-165
Praesto® Jetted Q35	1 L	PR00436-310
Praesto® Jetted Q35	5 L	PR00436-311
Praesto® Jetted Q35	10 L	PR00436-312

PRODUCT

COLUMN SIZE

ORDER NUMBER

Praesto® Jetted Q35 HT	1 ml Column	PR00436-275*
Praesto® Jetted Q35 HT	5 ml Column	PR00436-276*

** 1 ml HT columns packed with Ion Exchange resins available in packs of 5 only*

Contact information

Americas

150 Monument Road
Bala Cynwyd, PA
19004

T +1 800.343.1500
T +1 610.668.9090
F +1 484.384.2751

Americas@purolite.com

Europe

Llantrisant Business Park
Llantrisant
Wales, UK
CF72 8LF

T +44 1443 229334
F +44 1443 227073

Europe@purolite.com

Asia Pacific

Room 707, C Section
Huanglong Century Plaza
No.3 Hangda Road
Hangzhou, Zhejiang,
China 310007

T +86 571 876 31382
F +86 571 876 31385

AsiaPacific@purolite.com



Purolite®
Life Sciences

Algeria
Australia
Bahrain
Brazil
Canada
China
Czech Republic
France
Germany
India

Indonesia
Israel
Italy
Japan
Jordan
Kazakhstan
Korea
Malaysia
Mexico
Morocco

New Zealand
Poland
Romania
Russia
Singapore
Slovak Republic
South Africa
Spain
Taiwan
Tunisia

Turkey
UK
Ukraine
USA
Uzbekistan



Purolite Life Sciences brings Purolite's innovative thinking and distinguished history of resin technology expertise to the global Life Sciences marketplace.

Over three decades, Purolite has grown into the world's premier resin technology manufacturer and innovation leader, with production plants and advanced research labs across the globe.



100% focused on resin technology.



Global manufacturing at facilities
in the UK, Romania, China and USA.



De-risked long-term supply
through dual-sourcing.




25+ years of regulatory experience
from FDA inspected cGMP facility.



Over 35+ years of experience in solving
advanced R&D and purification challenges.





Purolite[®]
Life Sciences