

Purolite Ltd. Llantrisant Business Park Llantrisant South Wales CF72 8LF

T: +44 1443 229334 F: +44 1443 227073



**Purolite Life Sciences** 

www.purolitelifesciences.com lifesciences@purolite.com



A high-performance Protein A chromatography resin for increased productivity in downstream processing of clinical trial material such as monoclonal antibodies. It is designed to be a direct drop-in replacement for MabSelect™ SuRe LX\*.

# **About Praesto AP+**

Praesto AP+ (also referred to as 'Praesto AP Plus') is our latest-generation affinity chromatography resin to be produced using batch emulsification technology. It is an alkaline-stable Protein A affinity resin based on the first-generation product, Praesto AP, with enhanced performance characteristics over its predecessor.

Praesto AP+ has been purpose-designed for the production of clinical trial material and commercial manufacturling by significantly increasing the flow properties of Praesto AP, whilst maintaining high key performance properties such as dynamic binding capacity, CIP stability and packing reproducibility.

## **Base matrix**

The Praesto high-flow agarose base matrix provides a rigid, yet open pore structure. Compared to softer cross-linked agarose beads with similar porosity characteristics, Praesto resins can be operated at higher flow velocities at process scale, increasing productivity significantly.

## Market-leading performance

Praesto AP+ performance is benchmarked against MabSelect™ SuRe LX, from Cytiva Life Sciences. It therefore achieves very similar results for key performance parameters. This makes Praesto AP+ a high-performance resin which is also suitable as a direct drop-in replacement for MabSelect™ SuRe LX.

## Advantages:

- **High binding capacity** 50 g lgG per L of resin at 6 minutes residence time
- **Excellent flow properties** Max. flow velocity of 400 cm/h with pressure drop of less than 1.5 bar
- Alkaline stable At over 100 hours exposure to 0.1 M NaOH
- Self-pack or pre-packed columns Our column packing experts can work with your scientists
- Secure, validated supply at short lead times 500 L orders delivered in <30 days

\*MabSelect is a registered trademark of Cytiva Life Sciences

Figure 1: Dynamic binding capacity of Praesto AP+, Praesto AP, and MabSelect™ SuRe LX




Figure 2: Alkaline stability of Praesto AP+





Figure 3: Pressure/flow performance of Praesto AP+ packed into a BPG 300 column at 25 cm, 20 cm and 15 cm bed heights.



Typical Physical and Chemical Characteristics	
Polymer Structure	Highly cross-linked agarose
Appearance	Spherical beads
<b>Functional Group</b>	Alkaline-stable Protein A
Particle Size (d50v) μm	85
Dynamic Binding Capacity*	50 mg/ml (human lgG)
Maximum Flow Velocity (cm/h)**	400
pH Stability	2 - 12
Storage	20% ethanol (2 - 8°C)

- \* Determined at 10% breakthrough using polyclonal human immunoglobulin G at a residence time of 6 minutes.
- \*\* Determined using demineralized water at 20°C in a BPG 300 column (30 cm Inner Diameter) at 20 cm bed height with a pressure drop of less than 1.5 bar.

# **Ordering Information**

To place your order contact us via email or telephone and quote the product codes found in the table.

As part of our ongoing commitment to improving security of supply of bioprocess resins globally, we aim to deliver Praesto Protein A products in volumes smaller than 500 litres within 30 days of order.

# **Security of Supply**

Purolite® has manufacturing facilities at 5 strategic locations in the USA, Asia and Europe. Should Praesto production facilities in the UK be adversely affected to a level that would disrupt business operations, Purolite has designed manufacturing processes that facilitate rapid deployment at a different location. In this event, a rebuild plan would be implemented at our facility in Romania using an identical, modular design of the UK-based production facility.

We're ready to solve your process challenges. For further information on Purolite products and services, visit www.purolite.com or contact your nearest Technical Sales Office.

Bulk Resin	
VOLUME	PRODUCT CODE
10 ml	PR00330-163
25 ml	PR00330-166
100 ml	PR00330-164
500 ml	PR00330-165
1L	PR00330-310
5 L	PR00330-311
10 L	PR00330-312
10 L - 500 L	CONTACT US

Pre-Packed Columns	
COLUMN / VOLUME	PRODUCT CODE
HT Column (1 ml)	PR00330-275
HT Column (5 ml)	PR00330-276
HT Column (5 x 1 ml)	PR00330-575
HT Column (5 x 5 ml)	PR00330-576
MiniChrom (1 ml)	PR00330-175
MiniChrom (5 ml)	PR00330-176
RoboColumn (200 μl)	PR00330-174
RoboColumn (600 μl)	PR00330-279



Algeria
Australia
Bahrain
Brazil
Canada
China
Czech Republic
France
Germany

India Indonesia Israel Italy Japan Jordan Kazakhstan Korea Malaysia Mexico Morocco New Zealand Poland Romania Russia Singapore Slovak Republic South Africa Spain Taiwan Tunisia Turkey UK Ukraine USA Uzbekistan



## **Americas**

Purolite Corporation 2201 Renaissance Blvd. King of Prussia, PA 19406 T+1 800 343 1500 T+1 610 668 9090 F+1 800 260 1065 americas@purolite.com

## **EMEA**

Purolite Ltd.
Unit D
Llantrisant Business Park
Llantrisant, Wales, UK
CF72 8LF
T +44 1443 229334
F +44 1443 227073
emea@purolite.com

### **FSU**

Purolite Ltd.
Office 6-1
36 Lyusinovskaya Str.
Moscow, Russia
115093
T+7 495 363 5056
F+7 495 564 8121
fsu@purolite.com

## Asia Pacific

Purolite China Co. Ltd.
Room 707, C Section
Huanglong Century Plaza
No.3 Hangda Road
Hangzhou, Zhejiang, China 310007
T+86 571 876 31382
F+86 571 876 31385
asiapacific@purolite.com

Purolite, the leading manufacturer of quality ion exchange, catalyst, adsorbent and specialty high-performance resins, is the only company that focuses 100% of its resources on the development and production of resin technology.

We're ready to solve your process challenges. For further information on Purolite products and services, visit www.purolite.com or contact your nearest Technical Sales Office.



www.purolite.com

©2021 Purolite Corporation All rights reserved.

The statements, technical information and recommendations contained herein are believed to be accurate as of the date hereof. Since the conditions and methods of use of the product and of the information referred to herein are beyond our control, Purolite expressly disclaims any and all liability as to any results obtained or arising from any use of the product or reliance on such information; NO WARRANTY OF FITNESS FOR ANY PARTICULAR PURPOSE, WARRANTY OF MERCHANTABILITY OR ANY OTHER WARRANTY, EXPRESSED OR IMPLIED, IS MADE CONCERNING THE GOODS DESCRIBED OR THE INFORMATION PROVIDED HEREIN. The information provided herein relates only to the specific product designated and may not be applicable when such product is used in combination with other relaterals or in any process. Nothing contained herein constitutes a license to practice under any patent and it should not be construed as an inducement to infringe any patent and the user is advised to take appropriate steps to be sure that any proposed use of the product will not result in patent infringement.