PRODUCT DATA SHEET

Chromalite® MDEA/C

Diethylamino Methacrylate Macroporous Weak Anion Exchanger The Chromalite® M range of chromatographic resins are hydrophilic, macroporous, methacrylic resins for large-scale applications. Based on a rigid polymeric backbone that ensures excellent pressure-flow properties, they have been designed for high performance, stability and reliability under the conditions of modern industrial bioprocesses.

Chromalite® M resins for ion exchange chromatography include both strong and weak functional groups in both anion and cation exchange processes.

Chromalite® MDEA/C is a weak anion exchanger with optimal porosity and high mechanical stability.

PRINCIPAL APPLICATIONS

- Ion exchange chromatography
- Preparative chromatography
- Purification of peptides obtained by fermentation
- Dairy protein purification and fractionation
- Protein purification (up to 250 kDa)
- Vegetable Protein Fractionation

ADVANTAGES

- Hydrophillic
- · High mechanical stability
- * Shipped in wet form

TYPICAL PACKAGING

- 25ml
- 100ml
- 500ml
- 1L
- 5L
- 10L

TYPICAL PHYSICAL & CHEMICAL CHARACTERISTICS:

Particle Size - µm	100 - 300
Functional Group	Diethylamino Methacrylate
Pore Diameter	1000 Å
Working pressure (max.)	20 bar
Total Capacity (min.)	20 mg/ml
Optimal storage condition	2 - 20 °C
Supplied as	20% EtOH slurry



Americas

T +1 610 668 9090 F +1 610 668 8139 americas@purolite.com **EMEA**

T +44 1443 229334 F +44 1443 227073 europe@purolite.com Asia Pacific

T +86 571 876 31382 F +86 571 876 31385 asiapacific@purolite.com