

Chromalite® MDEA/F

**Diethylamino Methacrylate
Macroporous Weak Anion
Exchanger**

The Chromalite® M range of chromatographic resins are hydrophilic, macroporous, methacrylic resins for large-scale applications. Based on a rigid polymeric backbone that ensures excellent pressure-flow properties, they have been designed for high performance, stability and reliability under the conditions of modern industrial bioprocesses.

Chromalite® M resins for ion exchange chromatography include both strong and weak functional groups in both anion and cation exchange processes.

Chromalite® MDEA/F is a weak anion exchanger with optimal porosity and high mechanical stability.

This product is manufactured using [Jetting technology](#), our patented process for uniform particle size beads for superior performance.

PRINCIPAL APPLICATIONS

- Ion exchange chromatography
- Vegetable Protein Fractionation
- Protein purification (up to 250 kDa)
- Dairy protein purification and fractionation
- Purification of peptides obtained by fermentation
- Preparative chromatography

TYPICAL PACKAGING

- 25ml
- 100ml
- 500ml
- 1L
- 5L
- 10L

ADVANTAGES

- Hydrophilic
- High mechanical stability
- Uniform particle size

** Shipped in 20% Ethanol*

TYPICAL PHYSICAL & CHEMICAL CHARACTERISTICS:

Particle Size - µm	40 - 90
Functional Group	Diethylamino Methacrylate
Pore Diameter	1000 Å
Working pressure (max.)	20 bar
Total Capacity (min.)	20 mg/ml
Optimal storage condition	2 - 20 °C
Supplied as	20% EtOH slurry

