## PRODUCT DATA SHEET

# Chromalite® MEP/M

**Epoxy Methacrylate Macroporous Affinity Resin** 

The Chromalite® M range of chromatographic resins are hydrophilic, macroporous, methacrylic resins for large-scale applications. Based on a rigid polymeric backbone that ensures excellent pressure-flow properties, they have been designed for high performance, stability and reliability under the conditions of modern industrial bioprocesses.

Chromalite MEP/M is an affinity chromatography resin used in protein purification.

#### PRINCIPAL APPLICATIONS

- Preparative chromatography
- Ligand coupling for medical device application
- enzyme immobilization
- Purification of peptides obtained by fermentation
- Dairy protein purification and fractionation
- Protein purification (up to 250 kDa)
- Vegetable Protein Fractionation

### **ADVANTAGES**

- · High chemical stability
- Hydrophillic
- \* Shipped in wet form

#### **TYPICAL PACKAGING**

- 25ml
- 100ml
- 500ml
- 1L
- 5L
- 10L

# TYPICAL PHYSICAL & CHEMICAL CHARACTERISTICS:

Particle Size - μm	75 - 125
Functional Group	Ероху
Pore Diameter	1000 Å
Working pressure (max.)	20 bar
Optimal storage condition	2 - 20 °C
Total Capacity (min.)	n/a mg/ml
Supplied as	20% EtOH slurry



Americas

T +1 610 668 9090 F +1 610 668 8139 americas@purolite.com **EMEA** 

T +44 1443 229334 F +44 1443 227073 europe@purolite.com Asia Pacific

T +86 571 876 31382 F +86 571 876 31385 asiapacific@purolite.com